



# 17-Alpha-Hydroxyprogesterone-3-cmo [KLH] (DAGA-158K)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	17-Alpha-Hydroxyprogesterone-3-CMO, a progestogen which may reduce the risk of preterm delivery in pregnant women.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Keywords** 17-ALPHA-HYDROXYPROGESTERONE-3-CMO;