



Recombinant Chikungunya E2 Antigen (DAGA-162)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Chikungunya is an infection caused by the chikungunya virus. The disease features the sudden onset of fever two to four days after exposure. The fever usually lasts two to seven days, while accompanying joint pains typically last weeks or months but sometimes years. The mortality rate is a little less than 1 in 1,000; the elderly or those with underlying chronic medical problems are most likely to have severe complications.
----------------------------	---

Conjugate	Unconjugated
------------------	--------------

Applications	WB standard, antibody ELISA, immunogen, etc.
---------------------	--

Concentration	Batch dependent - please inquire should you have specific requirements.
----------------------	---

Size	1 mg, 10 mg
-------------	-------------

Preservative	None
---------------------	------

BACKGROUND

Keywords	CHIKV;Chikungunya Virus;Chikungunya E1;Chikungunya E2;CHIKV E1;CHIKV E2
-----------------	---