



## Trenbolone [KLH] (DAGA-178K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Trenbolone is a steroid used by veterinarians on livestock to increase muscle growth and appetite. To increase its effective half-life, trenbolone is not used in an unrefined form, but is rather administered as ester derivatives such as trenbolone acetate.
<b>Conjugate</b>	KLH
<b>Applications</b>	Immunogen
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

**Keywords** Trenbolone; 17 $\beta$ -Hydroxyestra-4; 9; 11-trien-3-one