



Dicloxacillin [KLH] (DAGA-187K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Dicloxacillin is a narrow-spectrum beta-lactam antibiotic of the penicillin class. It is used to treat infections caused by susceptible Gram-positive bacteria. It is active against beta-lactamase-producing organisms such as *Staphylococcus aureus*, which would otherwise be resistant to most penicillins. It is very similar to flucloxacillin and these two agents are considered interchangeable.

Nature	Synthetic
Expression System	N/A
Conjugate	KLH
Applications	Immunogen
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Dicloxacillin; Diclocil;