



Flumequine [KLH] (DAGA-193K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Flumequine is a first generation drug of the fluoroquinolone group of antibiotics. It has general use in the treatment of enteric infections in domestic species.
Nature	Synthetic
Expression System	N/A
Conjugate	KLH
Applications	Immunogen
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords 9-Fluoro-6; 7-dihydro-5-methyl-1-oxo-1H; 5H-benzo[ij]-quinolizine-2-carboxylic acid; Apurone;