



Sulfadoxine [BSA] (DAGA-202B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Sulfadoxine (also spelled sulphadoxine) is an ultra-long-lasting sulfonamide often used in combination with pyrimethamine to treat or prevent malaria. It is also used, usually in combination with other drugs, to treat or prevent various infections in livestock.
Nature	Synthetic
Expression System	N/A
Conjugate	BSA
Applications	LFIA
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	sulphadoxine; sulphormethoxine; wr4873; SULFADOXIN; SULFADOXINE; sulformethoxine; sulforthomidine; 6-(4-aminobenzenesulfonamide)-4; 5-dimethoxypyrimidine; 4-sulfanilamido-5; 6-dimethoxypyrimidine
-----------------	---