



Sulfonamides [KLH] (DAGA-211K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Sulfonamide or sulphonamide is the basis of several groups of drugs. The original antibacterial sulfonamides (sometimes called sulfa drugs or sulpha drugs) are synthetic antimicrobial agents that contain the sulfonamide group.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Sulfonamides; 1; 3-Difluoro-5-methylsulfonylbenzene; 3; 5-Difluoro-1-(methylsulfonyl)benzene; 2-Chloro-1-fluoro-4-methylsulfonylbenzene; 3-Chloro-4-fluoro-1-(methylsulfonyl)benzene; N-(4-Methoxybenzyl)-1; 3-propanesultam; N-(4-Methoxybenzyl)-N-methylbenzenes
-----------------	--