



## DENV type 2 NS1 (a.a. 777-1131) (DAGA-274)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	NS1 is essential for viral viability but its precise biological function is unknown. Antibodies raised in response to NS1 in viral infection can cross react with cell surface antigens on epithelial cells and platelets and this has been implicated in the d
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Suitable for use in ELISA and Lateral Flow. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Molecular Weight</b>	50 kDa
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Preservative</b>	None
<b>Storage</b>	Store at < -20°C.

### BACKGROUND

<b>Keywords</b>	Dengue Virus 1+2+3+4;Dengue Virus Type 1, 2, 3, 4;Group IV;Flaviviridae;Flavivirus;Dengue virus;Dengue Virus;DENV
-----------------	-------------------------------------------------------------------------------------------------------------------