



## Maduramicin [BSA] (DAGA-290B)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

**Antigen Description** Maduramicin is a monoglycoside polyether derived from the fungus *Actinomadura rubra*. It is classified as an ionophor and predominantly used as an anti-coccidial agent, which is mixed with poultry feed in the ratio of 6ppm.

**Conjugate** BSA

**Applications** LFIA

**Format** Liquid

**Concentration** Batch dependent - please inquire should you have specific requirements.

**Size** 1 mg

**Buffer** 0.01M pH7.4 PBS

**Preservative** None

**Storage** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

**Keywords** Maduramicin; Maduramicin A