



Native Dengue protein (Subtype 2) (DAGA-3010)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	There are four types of dengue viruses that cause dengue fever worldwide - Dengue 1, 2, 3 and 4. A person infected with one type of dengue will subsequently only be immune to that type. They will not be immune to other types of dengue and will, in fact, be at risk of developing severe symptoms if they contract another type of dengue.
Conjugate	Unconjugated
Applications	WB standard, Ab ELISA, immunogen, etc.
Format	Solution
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Supplied in PBS buffer, pH 7.4 with <0.1% CHAPS.
Storage	Aliquot and store at -70 centigrade.