



Recombinant Astrovirus Capsid protein [GST] (DAGA-3034)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Astrovirus Capsid Recombinant Product does contain a GST fusion partner.
Purity	Multi-step procedure including affinity chromatography.
Conjugate	GST
Applications	Suitable for use as a control in Lateral Flow Assays and as an immunogen. Suitable for spiking into stool samples. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Molecular Weight	60 kDa
Format	Purified
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	0.1 mg
Buffer	Phosphate Buffered Saline
Preservative	0.09% Sodium Azide
Storage	Store at 2-8°C for < 1 week. For long term storage, aliquot and store at -20°C to avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords Astrovirus; Astrovirus Capsid protein; Astrovirus Capsid