



Recombinant Norovirus GI.1 P Domain [His] (DAGA-3035)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Norovirus GI.1 P Domain Recombinant Product contains a fusion partner: His-tag.
Purity	Multi-step procedure including affinity chromatography.
Conjugate	His
Applications	Suitable for use as a control in Lateral Flow Assays and as an immunogen. Suitable for spiking into stool samples. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Molecular Weight	34 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	100 µg
Buffer	Tris Buffer
Preservative	0.09% Sodium Azide
Storage	Store at 2-8°C.