



# Recombinant EBV gp350/220 (RBD, a.a. 4-450) [His] (DAGA-3072)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	C-terminal 6xHis tagged gp350/220 (RBD/(a.a.4-450))(EBV) (a.a.4-450) (Genebank #: AHA36359). The receptor binding region of the gp350 protein.
<b>Purity</b>	>=95%
<b>Conjugate</b>	His
<b>Applications</b>	WB standard, antibody ELISA, immunogen, etc
<b>Format</b>	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	100 µg
<b>Storage</b>	Keep it at 4 centigrade if used within a month. For long term storage, split it into small aliquots and keep at -80 centigrade. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.