



Recombinant EBV gp350/220 (Ectodomain) [His] (DAGA-3073)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	C-terminal 6xHis tagged gp350/220 (RBD/(a.a.4-863))(EBV)(Genebank #: AHA36359).
Purity	>=95%
Conjugate	His
Applications	WB standard, antibody ELISA, immunogen, etc
Format	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	100 µg
Storage	Keep it at 4 centigrade if used within a month. For long term storage, split it into small aliquots and keep at -80 centigrade. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.