



Recombinant EBV p18 (DAGA-3074)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity>95%, based on SDS PAGEConjugateTBDApplicationsWB standard, antibody ELISA, immunogen, etc.FormatEach vial contains 100 μg of lyophilized protein in 50mM Tris-HCl pH 8.0, 60mM NaCl and 50% glycerol. Reconstitute the protein with 100 μl of Millipore water.ConcentrationBatch dependent - please inquire should you have specific requirements.Size100 μgStorageBefore reconstitution, stable for 1 year at -20 centigrade from the date of shipment. After reconstitution, stable for a month at 4 centigrade. Nonhazardous. No MSDS required.	Product Overview	The E.Coli derived recombinant mosaic protein contains the HHV-4 p18 regions, 1-119 amino acids.
Applications WB standard, antibody ELISA, immunogen, etc. Format Each vial contains 100 μg of lyophilized protein in 50mM Tris-HCl pH 8.0, 60mM NaCl and 50% glycerol. Reconstitute the protein with 100 μl of Millipore water. Concentration Batch dependent - please inquire should you have specific requirements. Size 100 μg Storage Before reconstitution, stable for 1 year at -20 centigrade from the date of shipment. After	Purity	>95% , based on SDS PAGE
Format Each vial contains 100 μg of lyophilized protein in 50mM Tris-HCl pH 8.0, 60mM NaCl and 50% glycerol. Reconstitute the protein with 100 μl of Millipore water. Concentration Batch dependent - please inquire should you have specific requirements. Size 100 μg Storage Before reconstitution, stable for 1 year at -20 centigrade from the date of shipment. After	Conjugate	TBD
glycerol. Reconstitute the protein with 100 µl of Millipore water. Concentration Batch dependent - please inquire should you have specific requirements. Size 100 µg Storage Before reconstitution, stable for 1 year at -20 centigrade from the date of shipment. After	Applications	WB standard, antibody ELISA, immunogen, etc.
Size 100 μg Storage Before reconstitution, stable for 1 year at -20 centigrade from the date of shipment. After	Format	Each vial contains 100 μg of lyophilized protein in 50mM Tris-HCl pH 8.0, 60mM NaCl and 50% glycerol. Reconstitute the protein with 100 μl of Millipore water.
Storage Before reconstitution, stable for 1 year at -20 centigrade from the date of shipment. After	Concentration	Batch dependent - please inquire should you have specific requirements.
	Size	100 μg
	Storage	

Email: info@creative-diagnostics.com

Tel: 1-631-624-4882 Fax: 1-631-938-8221