



# Recombinant Enterovirus B Echovirus 18 [GST] (DAGA-3076)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	>85% as determined by SDS-PAGE
Conjugate	GST
Applications	IM, LF, ST
Molecular Weight	53 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	0.1 mg
Buffer	Tris buffer with salts, pH 8.0.
Preservative	0.09% Sodium Azide
Storage	Store at 2-8°C for 1 week. For long term storage, aliquot and store at -20°C to avoid multiple freeze/thaw cycles.

## BACKGROUND

Introduction	Echovirus is a polyphyletic group of viruses associated with enteric disease in humans. The name is derived from "enteric cytopathic human orphan virus". These viruses were originally not associated with disease, but many have since been identified as disease-causing agents. The term "echovirus" was used in the scientific names of numerous species, but all echoviruses are now recognized as strains of various species, most of which are in the family Picornaviridae.
--------------	--

**Keywords**

Enterovirus; EV; Enterovirus B Echo 18; EV B Echo 18

---