



Recombinant Enterovirus C Coxsackievirus A24 [GST] (DAGA-3077)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	>85% as determined by SDS-PAGE
Conjugate	GST
Applications	IM, LAT, ST
Molecular Weight	61 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	0.1 mg
Buffer	Tris buffer with salts, pH 8.0.
Preservative	0.09% Sodium Azide
Storage	Store at 2-8°C for 1 week. For long term storage, aliquot and store at -20°C to avoid multiple freeze/thaw cycles.

BACKGROUND

Introduction	Coxsackievirus is a virus that belongs to a family of nonenveloped, linear, positive-sense single-stranded RNA viruses, Picornaviridae and the genus Enterovirus, which also includes poliovirus and echovirus. Enteroviruses are among the most common and important human pathogens, and ordinarily its members are transmitted by the fecal-oral route.
Keywords	Enterovirus; EV; Enterovirus C Cox A24; EV C Cox A24