



# Recombinant Dog allergen, can f5. [His] (DAGA-3086)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Nature</b>	Recombinant
<b>Expression System</b>	Pichia pastoris
<b>Purity</b>	> 95% pure (SDS-PAGE).
<b>Conjugate</b>	His
<b>Applications</b>	It's being tested in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Purified
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	0.5 mg
<b>Buffer</b>	20 mM Phosphate Buffer, pH 8, 5 mM EDTA and 0.15 M Sodium Chloride.
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C to -80°C. Aliquot and store at -80°C to avoid multiple freeze/thaw cycles.