



## Toltrazuril [BSA] (DAGA-337B)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Antigen Description	Toltrazuril is a coccidiostat.
Nature	Synthetic
Expression System	N/A
Conjugate	BSA
Applications	LFIA
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

Keywords	Toltrazuril; Baycox;
----------	----------------------