



Methotrexate [KLH] (DAGA-351K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Methotrexate is in the class of medications called antimetabolites. The drug is used to treat severe psoriasis, rheumatoid arthritis and certain types of cancer. It works by slowing the proliferation of cells. Specifically, it inhibits the enzyme dihydrofolate reductase that is necessary for the synthesis of folic acid, which is required for the synthesis of DNA, RNA, thymidylates and proteins.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	L-(+)-AMETHOPTERIN; Methotrexate; MTX; ANTIFOLAN; NSC 740; 4-AMINO-N 10-METHYLPTEROYL-L-GLUTAMIC ACID; 4-AMINO-10-METHYLPTEROYLGLUTAMIC ACID;
-----------------	---