



Recombinant Protein *Akkermansia muciniphila* (strain ATCC BAA-835 / Muc) atpF, Yeast (DAGA-3633)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The following tags are available. N-terminal His-tagged Tag-Free
Antigen Description	Recombinant 6xHis tagged <i>Akkermansia muciniphila</i> (strain ATCC BAA-835 / Muc) atpF (UniProtKB No. B2UNX8) was expressed.
Species	<i>Akkermansia muciniphila</i>
Purity	>85% (SDS-PAGE)
Conjugate	His
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	Tris-based buffer, 50% glycerol
Storage	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.