



# Adenovirus type 5 hexon protein (>99%) (DAGA-450)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Adenoviruses are DNA viruses generally widespread in nature that are frequently the cause of acute upper respiratory tract infections (i.e. common colds). Forty-seven known serotypes have been isolated since they were first discovered in 1953 with 4 types
<b>Purity</b>	>99% by SDS-PAGE
<b>Conjugate</b>	Unconjugated
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	100 µg
<b>Buffer</b>	10mM Bis Tris Propane, 375 mM NaCl
<b>Preservative</b>	None
<b>Storage</b>	Store at -80°C

## BACKGROUND

**Keywords** Adenovirus type 5 hexon antigen; Adenovirus type 5 protein; Adenovirus hexon protein