



Recombinant Parvovirus B19 VP2 VLP (a.a. 228-781) (>90%) (DAGA-519)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Parvoviruses cause infections in a wide variety of birds and mammals. Parvovirus, commonly called parvo, is a genus of the Parvoviridae family linear, non-segmented single stranded DNA viruses with an average genome size of 9 kbp.
Purity	>90% by SDS-PAGE
Conjugate	TBD
Applications	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Molecular Weight	65 kD
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Buffer	Phosphate Buffered Saline, pH 7.4, 10% Glycerol
Preservative	None
Storage	Store at -65°C. Aliquot to avoid repeated freeze-thaw cycles

BACKGROUND

Keywords Parvovirus protein;Parvo protein;Parovirus;parvovirus protein