



# Recombinant parvovirus B19 VP1 (96%) [His] (DAGA-523)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Parvoviruses cause infections in a wide variety of birds and mammals. Parvovirus, commonly called parvo, is a genus of the Parvoviridae family linear, non-segmented single stranded DNA viruses with an average genome size of 13 kbp.
<b>Purity</b>	>96% by SDS-PAGE
<b>Conjugate</b>	His
<b>Applications</b>	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Molecular Weight</b>	55 kD
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Buffer</b>	50mM Tris/HCl pH 8.0, 500mM Sodium Chloride and 8M Urea.
<b>Preservative</b>	None
<b>Storage</b>	Store at <-70°C.

## BACKGROUND

**Keywords** Parvovirus protein;Parvo protein;Parovirus;parvovirus protein