



# Recombinant Parvovirus B19 Pep A Protein (aa 1-230) (DAGA-527)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Parvoviruses have a single kind of coat subunit, but a few of the 60 copies of this protein have an extra,N-terminal domain, which is sequestered within the virion.
<b>Nature</b>	Recombinant
<b>Expression System</b>	E. coli
<b>Conjugate</b>	TBD
<b>Molecular Weight</b>	26 kDa
<b>Procedure</b>	None
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	1 mg
<b>Buffer</b>	50 mM Tris/HCl pH 8.0, 500 mM NaCl and 8 M urea
<b>Preservative</b>	None
<b>Storage</b>	-65°C or lower

## BACKGROUND

<b>Keywords</b>	Parvoviridae;Parvovirus B19;Parvovirus B19;VP1;Parvovirus (B19);B19 virus;Parvovirus B19 VP1
-----------------	--