



Recombinant Parvovirus B19 VP2 Protein (aa 233-781) (DAGA-531)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Parvoviruses have a single kind of coat subunit, but a few of the 60 copies of this protein have an extra,N-terminal domain, which is sequestered within the virion.
Conjugate	TBD
Molecular Weight	63 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	50 mM Tris/HCl pH 8.0, 500 mM NaCl and 8 M urea.
Preservative	None
Storage	-65°C or lower

BACKGROUND

Keywords	Parvoviridae;Parvovirus B19;Parvovirus B19;VP2;Parvovirus (B19);B19 virus;Parvovirus B19 VP2
----------	--