



# Recombinant Parvovirus B19 NS1 Protein (DAGA-579)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	The parvoviruses have small ( 25 nm diameter) and structurally stable capsids that bind host cell receptors, small ssDNA genomes with few genes that control their interactions with their hosts, and they replicate only in dividing cells or in the presence of a helper virus.
<b>Conjugate</b>	TBD
<b>Molecular Weight</b>	78 kDa
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	1 mg
<b>Buffer</b>	50 mM Tris/HCl pH 8.0, 500 mM NaCl and 8 M urea.
<b>Preservative</b>	None
<b>Storage</b>	-65°C or lower

## BACKGROUND

<b>Keywords</b>	Parvoviridae;Parvovirus B19;Parvovirus B19;NS1;Parvovirus (B19);B19 virus;Parvovirus B19 NS1
-----------------	--