



# Recombinant SARS Nucleoprotein (a.a. 1-49, 192-220) (DAGA-606)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	E.coli derived recombinant. The protein contains Nucleocapsid protein immunodominant regions
<b>Purity</b>	>95% by SDS-PAGE
<b>Conjugate</b>	TBD
<b>Applications</b>	Antigen in ELISA and Western blots, excellent antigen for detection of SARS with minimal specificity problems.
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	1 mg
<b>Buffer</b>	50 mM Tris-HCl, 60 mM Na Cl, 50% glycerol.
<b>Preservative</b>	None
<b>Storage</b>	Upon arrival, store at -20°C

## BACKGROUND

<b>Keywords</b>	Coronavirus Core protein;SARS;Core protein;Core protein SARS
-----------------	--