



HHV-4 p18 (DAGA-612)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	E.coli derived recombinant. The mosaic protein contains the HHV-4 p18 regions
Purity	>95% by SDS-PAGE
Conjugate	TBD
Applications	Antigen in ELISA and Western blots, excellent antigen for detection of HHV-4 (EBV)- with minimal specificity problems.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	25mM Tris-HCl, 1,5 M urea, 50% glycerol
Preservative	None
Storage	Protein is shipped at ambient temperature. Upon arrival, Store at -20°C

BACKGROUND

Keywords	Herpesviridae;p18;Lymphocryptovirus p18;Human herpesvirus 4 (HHV-4);p18
----------	---