



Dog albumin reference serum (DAGA-622)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Albumin refers generally to any protein that is water soluble, which is moderately soluble in concentrated salt solutions, and experiences heat denaturation. They are commonly found in blood plasma, and are unique to other plasma proteins in that they are not glycosylated. Substances containing albumin, such as egg white, are called albuminoids.
Specificity	Canine, Dog
Nature	Native
Expression System	Dog
Conjugate	Unconjugated
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords Dog;Albumin;Serum