



Dog immunoglobulin reference serum (DAGA-628)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Immunoglobulin (Ig), There are four IgG subclasses (IgG1, 2, 3 and 4) in humans, named in order of their abundance in serum (IgG1 being the most abundant). IgG2 is the only IgG subclass which passes through the placenta at a level generally lower than that found in the mother. A deficiency of IgG2 indicates a poor antibody response to bacterial polysaccharides and can lead to increased susceptibility to infections caused by encapsulated bacteria
Specificity	Canine, Dog
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords	Dog; Immunoglobulin; Serum
-----------------	----------------------------