



Dog retinol binding protein reference serum (DAGA-631)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Retinol-binding proteins (RBP) are a family of proteins with diverse functions. They are carrier proteins that bind retinol. Assessment of retinol-binding protein is used to determine visceral protein mass in health-related nutritional studies.
Specificity	Canine, Dog
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords	Dog;Retinol Binding Protein;Serum
-----------------	-----------------------------------