



Mouse alpha 1-antitrypsin reference serum (DAGA-656)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Alpha-1 Antitrypsin is a protease inhibitor belonging to the serpin superfamily. It is known as serum trypsin inhibitor. Alpha-1 antitrypsin is also referred to as alpha-1 proteinase inhibitor (A1PI) because it inhibits a wide variety of proteases. It protects tissues from enzymes of inflammatory cells, especially neutrophil elastase.
Specificity	Mouse
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords Mouse;Alpha 1-Antitrypsin;Serum