



Mouse retinol binding protein reference serum (DAGA-670)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|--|
| Antigen Description | Retinol-binding proteins (RBP) are a family of proteins with diverse functions. They are carrier proteins that bind retinol. Assessment of retinol-binding protein is used to determine visceral protein mass in health-related nutritional studies. |
| Specificity | Mouse |
| Conjugate | Unconjugated |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements |
| Size | 1 ml |
| Preservative | 0.1% Sodium Azide |
| Storage | Frozen -20°C |

BACKGROUND

Keywords Mouse;Retinol Binding Protein;Serum