



Rabbit haptoglobin reference serum (DAGA-680)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Haptoglobin (Hp) is an acute phase protein that binds the free hemoglobin (Hb), thus preventing iron loss and renal damage. Hp also has antioxidative and immunomodulatory properties. The haptoglobin assay is used to screen for and monitor intravascular hemolytic anemia. In intravascular hemolysis, free hemoglobin will be released into circulation and hence haptoglobin will bind the hemoglobin.

Specificity	Rabbit
Nature	Native
Expression System	Rabbit
Conjugate	Unconjugated
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords Rabbit;Haptoglobin;Serum