



# Rat alpha 1-acid glycoprotein reference serum (DAGA-684)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

**Antigen Description** Alpha-1-acid glycoprotein ( $\alpha$ 1AGp, AGP or AAG) is an acute phase (acute phase protein) plasma alpha-globulin glycoprotein and is modulated by two polymorphic genes. It is synthesized primarily in hepatocytes and has a normal plasma concentration between 0.6-1.2 mg/mL (1-3% plasma protein). Plasma levels are affected by pregnancy, burns, certain drugs, and certain diseases, particularly HIV.

<b>Specificity</b>	Rat
<b>Nature</b>	Native
<b>Expression System</b>	Rat
<b>Conjugate</b>	Unconjugated
<b>Procedure</b>	None
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	1 ml
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Frozen -20°C

## BACKGROUND

**Keywords** Rat;Alpha 1-Acid Glycoprotein;Serum