



# Rat beta 2-microglobulin reference serum (DAGA-686)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	$\beta$ 2 microglobulin also known as B2M is a component of MHC class I molecules, which are present on all nucleated cells (excludes red blood cells). In humans, the $\beta$ 2 microglobulin protein is encoded by the B2M gene.
<b>Specificity</b>	Rat
<b>Nature</b>	Native
<b>Expression System</b>	Rat
<b>Conjugate</b>	Unconjugated
<b>Procedure</b>	None
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	1 ml
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Frozen -20°C

## BACKGROUND

**Keywords** Rat;Beta 2-Microglobulin;Serum