



Rat hemopexin reference serum (DAGA-694)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Hemopexin (or haemopexin; HPX), also known as beta-1B-glycoprotein is a protein that in humans is encoded by the HPX gene and belongs to hemopexin family of proteins. Hemoglobin and its scavenger protein hemopexin (Hx) associate with HDL and influence the inflammatory properties of HDL. In addition it can also be said that HDL from Hx-null mice is proinflammatory. Moreover, Hemopexin deficiency is associated with various other inflammatory diseases such as septic shock and experimental autoimmune encephalomyelitis.

Specificity Rat

Conjugate Unconjugated

Format Liquid

Concentration Batch dependent - please inquire should you have specific requirements

Size 1 ml

Preservative 0.1% Sodium Azide

Storage Frozen -20°C

BACKGROUND

Keywords Rat;Hemopexin;Serum