



Rat SAP reference serum (DAGA-702)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|---|
| Antigen Description | The serum amyloid P component (SAP) is the identical serum form of amyloid P component (AP), a 25kDa pentameric protein first identified as the pentagonal constituent of in vivo pathological deposits called "amyloid". |
| Specificity | Rat |
| Conjugate | Unconjugated |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements |
| Size | 1 ml |
| Preservative | 0.1% Sodium Azide |
| Storage | Frozen -20°C |

BACKGROUND

| | |
|-----------------|---------------|
| Keywords | Rat;SAP;Serum |
|-----------------|---------------|