



Rat SAP reference serum (DAGA-702)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The serum amyloid P component (SAP) is the identical serum form of amyloid P component (AP), a 25kDa pentameric protein first identified as the pentagonal constituent of in vivo pathological deposits called "amyloid".
Specificity	Rat
Nature	Native
Expression System	Rat
Conjugate	Unconjugated
Procedure	None
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords	Rat;SAP;Serum
-----------------	---------------