



## V. cholerae Cholera Toxin B [Biotin] (DAGB151)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Species	V. cholerae
Conjugate	Biotin
Applications	ELISA
Reconstitution	When reconstituted with 0.4 ml of water the protein is in 0.01M sodium phosphate, pH7.5.
Endotoxin	The endotoxin content, determined using a kinetic chromagenic LAL assay, is <10 EU/mg.
Format	Lyophilized power, sealed under vacuum
Size	0.2 mg
Preservative	None
Storage	Store at 2-8°C prior to and after reconstitution.