



Eastern Equine Encephalitis Virus E3E2 Protein [His] (DAGB175)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Nature	Recombinant
Expression System	Insect cells
Species	Eastern Equine Encephalitis Virus
Conjugate	His
Applications	WB
Procedure	None
Format	Liquid
Concentration	0.901 mg/mL
Size	100 µg
Buffer	PBS (supplemented with glycerol, arginine and glutamic acid)
Preservative	None
Storage	2-3 weeks at -20°C, -80°C long term