



# Ebola virus Glycoprotein (minus the Transmembrane Region) [HA] (DAGB181)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Ebola virus
Conjugate	HA
Applications	ELISA, WB
Molecular Weight	~80-120 kDa including the HA-tag, without glycosylation
Format	Liquid
Concentration	1.008 mg/mL
Size	100 µg
Buffer	PBS
Preservative	None
Storage	2-3 weeks at -20°C, -80°C long term