



Robenidine [HRP] (DAGB311)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Each conjugate comprises antigen covalently bound to horseradish peroxidase and is suitable as a tracer in immunoassay development.
Species	N/A
Conjugate	HRP
Applications	IA
Format	The conjugate is supplied as a concentrate. Dilute as required and use working strength conjugate immediately after dilution
Size	0.5 ml
Preservative	None
Storage	2 - 8°C for up to 3 months / -20°C for long term storage

BACKGROUND

Introduction	Robenidine is a coccidiostat. Robenidine is an antibiotic used for the control of coccidiosis, a debilitating protozoal infection in poultry. Although there are alternative antibiotics available, robenidine is important in the management of antibiotic resistance as farmers rotate the use of robenidine with other antibiotics to try to preserve the effectiveness of these products in fighting infections.
---------------------	--