



## Fumonisin [HRP] (DAGB357)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Specificity</b>	Each conjugate comprises antigen covalently bound to horseradish peroxidase and is suitable as a tracer in immunoassay development.
<b>Species</b>	N/A
<b>Conjugate</b>	HRP
<b>Applications</b>	IA
<b>Format</b>	The conjugate is supplied as a concentrate. Dilute as required and use working strength conjugate immediately after dilution
<b>Size</b>	1 mg
<b>Preservative</b>	None
<b>Storage</b>	2 - 8°C for up to 3 months / -20°C for long term storage

### BACKGROUND

<b>Introduction</b>	The fumonisins are a group of mycotoxins derived from Fusarium, Liseola section. They have strong structural similarity to sphinganine, the backbone precursor of sphingolipids
---------------------	---