



Chicken Sterile Filtered Serum (DAGC282)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The raw materials for these products were collected from chickens in the USA in abattoirs registered with the USDA. The birds received ante and post-mortem inspections under a veterinarian's supervision and were apparently free from infectious and contagious diseases. At no time during collection or processing was the material commingled with any other material of animal origin.
Nature	Native
Expression System	Chicken serum
Species	Chicken
Purity	0.2 micron (Absolute)
Conjugate	Unconjugated
Applications	ELISA
Procedure	None
Format	Liquid
Concentration	Hemoglobin: <110 mg/dl Protein:>2.0 g/dl
Size	100 ml, 500 ml
Buffer	pH: 6.8 - 8.2
Preservative	None
Storage	This serum is stable for ≥ 3 years when stored $-20^{\circ}\text{C} \pm 5^{\circ}\text{C}$. Avoid repeated freezing and thawing.

BACKGROUND

Keywords

Chicken Serum; Serum
