



Equine Sterile Filtered Serum (DAGC286)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The serum was produced from animals in the USA. This serum was collected in USDA approved collection facilities. The animals received inspections under a veterinarian's supervision and were apparently free from infectious and contagious diseases. At no time during collection or processing was the material commingled with any other material of animal origin.
Nature	Native
Expression System	Donor Equine Serum
Species	Equine
Purity	0.2 micron (Absolute)
Conjugate	Unconjugated
Applications	ELISA
Procedure	None
Format	Liquid
Concentration	Hemoglobin: <20 mg/dl Protein: ≥5.0 g/dl
Size	100 ml, 500 ml
Buffer	pH: 6.8 - 8.2
Preservative	None
Storage	This serum is stable for ≥ 5 years when stored -20°C ± 5°C. Avoid repeated freezing and thawing.

BACKGROUND

Keywords

Equine Serum; Serum
