



# Fetal Bovine Sterile Filtered Serum (DAGC287)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

### Product Overview

The raw materials for these products were collected from animals in the USA in abattoirs registered with the USDA. The animals received ante and post-mortem inspections under a veterinarian's supervision and were apparently free from infectious and contagious diseases. At no time during collection or processing was the material commingled with any other material of animal origin.

Species	Bovine
Conjugate	Unconjugated
Applications	ELISA
Endotoxin	$\leq 10$ EU/ml
Format	Liquid
Concentration	Hemoglobin: $\leq 20$ mg/dl Protein: $\geq 3$ g/dl
Size	10 ml, 50 ml, 100 ml, 500 ml
Buffer	pH: 7.0-8.0
Preservative	None
Storage	For long term storage, store frozen at $-20^{\circ}\text{C} \pm 5^{\circ}\text{C}$ . Avoid repeated freezing and thawing.

## BACKGROUND

**Keywords** Fetal Bovine Serum; Serum; Bovine Serum