



Mouse Sterile Filtered Serum (DAGC289)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The raw material used was collected from Non-Swiss Albino mice. The material for this product was processed in the USA from the blood of healthy mice that had not been inoculated with or exposed to any infectious agents. At no time during collection or processing was the material commingled with any other material of animal origin
Nature	Native
Expression System	Mice, sexually mature, mixed sex, Non Swiss Albino
Species	Mouse
Purity	0.2 micron (Absolute)
Conjugate	Unconjugated
Applications	ELISA
Procedure	None
Format	Amber to slightly reddish color
Concentration	Hemoglobin: <35 mg Protein: ≥4.5 g/dl
Size	50 ml, 100 ml, 500 ml
Buffer	pH: 6.8 - 8.2
Preservative	None
Storage	This serum is stable for ≥ 3 years when stored -20°C ± 5°C. Avoid repeated freezing and thawing.

BACKGROUND

Keywords

Mouse Serum; Serum
