



# Rhesus Monkey Sterile Filtered Serum (DAGC293)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	The material for this product was processed in the USA from the blood of healthy Rhesus monkeys that had not been inoculated with or exposed to any infectious agents. At no time during collection or processing was the material commingled with any other material of animal origin. This product is not to be used for direct or indirect exposure or inoculation into laboratory and domestic animals.
<b>Nature</b>	Native
<b>Expression System</b>	Rhesus Monkey serum, mixed sex, mixed breed
<b>Species</b>	Rhesus Monkey
<b>Purity</b>	0.2 micron (Absolute)
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA
<b>Procedure</b>	None
<b>Format</b>	Liquid
<b>Concentration</b>	Hemoglobin: <30 mg/dl Protein:>4.0 g/dl
<b>Size</b>	50 ml, 100 ml, 500 ml
<b>Buffer</b>	pH: 6.8 - 8.2
<b>Preservative</b>	None

**Storage**

This serum is stable for  $\geq 3$  years when stored  $-20^{\circ}\text{C} \pm 5^{\circ}\text{C}$ . Avoid repeated freezing and thawing.

---

## BACKGROUND

**Keywords**

Rhesus Monkey Serum; Serum; Monkey Serum

---